


Product Datasheet

FOXA1 Antibody (orb1252290)


Description	FOXA1 Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	A human partial recombinant protein corresponding to amino acids 372-472 was used as the immunogen for this FOXA1 antibody.
Target	FOXA1
Preservatives	PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide
Form/Appearance	Liquid
Concentration	0.2 mg/mL
Storage	Aliquot and Store at 2-8°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	Flow Cytometry: 0.5-1ug/10 ⁶ cellsIF: 1-2 ug/mlWB: 0.5-1 ug/mlIHC (FFPE): 1-2 ug/ml for 30 min at RTOptimal dilution of the FOXA1 antibody should be determined by the researcher.
Isotype	IgG1, kappa
Clonality	Monoclonal
Uniprot ID	P55317
Dilution Range	Flow Cytometry: 0.5-1ug/10 ⁶ cellsIF: 1-2 ug/mlWB: 0.5-1 ug/mlIHC (FFPE): 1-2 ug/ml for 30 min at RTOptimal dilution of the FOXA1 antibody should be determined by the researcher.
Expiration Date	12 months from date of receipt.



IHC testing of
FFPE human
prostate
carci...



IHC testing of
FFPE human
prostate with
...



Western blot
testing of
human 1)
MCF7 (b...