


Product Datasheet

ATP2B1 antibody (orb125230)

Description	Goat Polyclonal to ATP2B1. NP_001001323.1; NP_001673.2, PMID: 21126504 and GeneID: 490.
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	ATP2B1 (aa312-327)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 150kDa band observed in lysates of cell line HeLa (calculated MW of 135kDa according to NP_001673.2). Recommended concentration: 0.03-0.1µg/ml.
Protein Sequence	KQDGA IENRNKAKAQD
Clonality	Polyclonal
MW	130; 135
Entrez	490
Dilution Range	ELISA: 1:64000, WB: 0.03-0.1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa lysate using ATP2B1 antibody (orb125230).