

Product Datasheet

DSG1 Antibody (orb1252309)

Description

DSG1 Antibody

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IF, IP, WB

Applications

Immunogen

A recombinant protein corresponding to the intracellular domain of human Desmoglein 1 was used as the immunogen for the Desmoglein antibody.

Target

DSG1

Preservatives

PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide

Form/Appearance

Liquid

Concentration

0.2 mg/mL

Storage

Aliquot and Store at 2-8°C. Avoid freez-thaw cycles.

Note

For research use only

Application notes

Flow Cytometry: 0.5-1 ug/million cells in 0.1ml
Immunofluorescence: 1-2 ug/ml
WB: 0.5-1 ug/ml
Immunoprecipitation: suitable
Optimal dilution of the Desmoglein antibody should be determined by the researcher.

Isotype

IgG1, kappa

Clonality

Monoclonal

Uniprot ID


[Q02413](#)

Dilution Range

Flow Cytometry: 0.5-1 ug/million cells in 0.1ml
Immunofluorescence: 1-2 ug/ml
WB: 0.5-1 ug/ml
Immunoprecipitation: suitable
Optimal dilution of the Desmoglein antibody should be determined by the researcher.

Expiration Date

12 months from date of receipt.



SDS-PAGE
Analysis of
Purified, BSA-
Free ...