


Product Datasheet

ZNF652 antibody (orb125233)

Description	Goat Polyclonal to ZNF652. NP_055712.1, PMID: 20204290 and GeneID: 22834.
Species/Host	Goat
Reactivity	Bovine, Human, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	ZNF652 (aa594-606)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 90kDa band observed in Human Heart, Kidney and Liver lysates (calculated MW of 69.7kDa according to NP_055712.1). The observed molecular weight corresponds to earlier findings in literature with different antibodies (Kumar et al, Mol Cancer Res. 20. IHC: Paraffin embedded Human Small Intestine. Recommended concentration: 5µg/ml.
Protein Sequence	RHLAEKNSSAQHH
Clonality	Polyclonal
MW	69.7
Entrez	22834
Dilution Range	ELISA: 1:64000, WB: 0.01-0.03 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Heart lys...