

Product Datasheet

Zfp652 antibody (orb125234)

Description	Goat Polyclonal to Zfp652. NP_963903.2, PMID: 20204290.
Species/Host	Goat
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA
Target	Zfp652 (mouse)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Preliminary experiments gave bands at approx 45kDa and 18kDa in Mouse Heart lysates after 1µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the bands we observe given the calculated size of 69.4kDa ac
Protein Sequence	RHLAEKNSAAPAQHH
Clonality	Polyclonal
MW	69.4
Dilution Range	ELISA: 1:16000, WB: 1 µg/ml
Expiration Date	12 months from date of receipt.