

Product Datasheet

Lgals1 Antibody (orb1252358)

Description

Lgals1 Antibody

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, FC, IF, IHC-P, WB

Applications

Immunogen

A human partial recombinant protein corresponding to amino acids 12-108 was used as the immunogen for the Galectin 1 antibody.

Target

Lgals1

Preservatives

PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide

Form/Appearance

Liquid

Concentration

0.2 mg/mL

Storage

Aliquot and Store at 2-8°C. Avoid freeze-thaw cycles.

Note

For research use only

Application notes

ELISA: 2-4 ug/ml; order the BSA-free formulation for coating
IF: 1-2 ug/ml
Flow Cytometry: 0.5-1 ug/million cells
in 0.1ml
WB: 0.5-1 ug/ml
IHC (FFPE): 0.5-1 ug/ml for 30 min at RT
titrating of the Galectin 1 antibody may be required for optimal performance.

Isotype

IgG1, lambda

Clonality

Monoclonal

Uniprot ID


[P09382](#)

Dilution Range


ELISA: 2-4 ug/ml; order the BSA-free formulation for coating
IF: 1-2 ug/ml
Flow Cytometry: 0.5-1 ug/million cells
in 0.1ml
WB: 0.5-1 ug/ml
IHC (FFPE): 0.5-1 ug/ml for 30 min at RT
titrating of the Galectin 1 antibody may be required for optimal performance.

Expiration Date

12 months from date of receipt.



IHC testing of
FFPE human
prostate
carci...



Western blot
testing of
human 1)
HeLa, 2...