



Product Datasheet

IFNG Antibody (orb1252361)


Description	IFNG Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Human recombinant protein was used as the immunogen for the IFN gamma antibody.
Target	IFNG
Preservatives	PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide
Form/Appearance	Liquid
Concentration	0.2 mg/mL
Storage	Aliquot and Store at 2-8°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	Flow Cytometry: 0.5-1 ug/million cells in 0.1ml Immunofluorescence: 1-2 ug/ml IHC (FFPE): 0.5-1 ug/ml for 30 min at RT Optimal dilution of the IFN gamma antibody should be determined by the researcher.
Isotype	IgG2b, kappa
Clonality	Monoclonal
Uniprot ID	P01579
Dilution Range	Flow Cytometry: 0.5-1 ug/million cells in 0.1ml Immunofluorescence: 1-2 ug/ml IHC (FFPE): 0.5-1 ug/ml for 30 min at RT Optimal dilution of the IFN gamma antibody should be determined by the researcher.
Expiration Date	12 months from date of receipt.



IHC testing of FFPE human renal cell carcinoma tissue.



IHC testing of FFPE human renal cell carcinoma tissue.



IHC testing of FFPE human testis tissue with IFN gamma antibody.