

Product Datasheet

CDH1 Antibody (orb1252364)

Description

CDH1 Antibody

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

FC, IF, IHC-P, WB

Immunogen

A full length recombinant human protein was used as the immunogen for the E-Cadherin antibody.

Target

CDH1

Preservatives

PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide

Form/Appearance

Liquid

Concentration

0.2 mg/mL

Storage

Aliquot and Store at 2-8°C. Avoid freeze-thaw cycles.

Note

For research use only

Application notes

Flow Cytometry: 0.5-1 ug/million cells in 0.1mMImmunofluorescence: 0.5-1 ug/mlWestern blot: 0.5-1 ug/mlImmunohistochemistry (FFPE): 1-2 ug/ml for 30 min at RT (1)Optimal dilution of the E-Cadherin antibody should be determined by the researcher.1. Staining of formalin-fixed tissues requires boiling tissue sections in 10mM Citrate buffer, pH 6.0, for 10-20 min followed by cooling at RT for 20 min.

Isotype

IgG1, kappa

Clonality

Monoclonal

Uniprot ID


P12830

Dilution Range


Flow Cytometry: 0.5-1 ug/million cells in 0.1mMImmunofluorescence: 0.5-1 ug/mlWestern blot: 0.5-1 ug/mlImmunohistochemistry (FFPE): 1-2 ug/ml for 30 min at RT (1)Optimal dilution of the E-Cadherin antibody should be determined by the researcher.1. Staining of formalin-fixed tissues requires boiling tissue sections in 10mM Citrate buffer, pH 6.0, for 10-20 min followed by cooling at RT for 20 min.

Expiration Date

12 months from date of receipt.



Formalin-fixed, paraffin-embedded human ...



Formalin-fixed, paraffin-embedded human ...



Formalin-fixed, paraffin-embedded human ...