


Product Datasheet

PECAM1 Antibody (orb1252365)

Description	PECAM1 Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Human full length recombinant protein was used as the immunogen for this CD31 antibody.
Target	PECAM1
Preservatives	PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide
Form/Appearance	Liquid
Concentration	0.2 mg/mL
Storage	Aliquot and Store at 2-8°C. Avoid freez-thaw cycles.
Note	For research use only
Application notes	Flow Cytometry: 0.5-1ug/10 ⁶ cellsIF: 1-2 ug/mlThe concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the CD31 antibody to be titrated up or down for optimal performance.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P16284
Dilution Range	Flow Cytometry: 0.5-1ug/10 ⁶ cellsIF: 1-2 ug/mlThe concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the CD31 antibody to be titrated up or down for optimal performance.
Expiration Date	12 months from date of receipt.



SDS-PAGE
Analysis of
Purified,
BSA-Free ...