

Product Datasheet

Macrophage / Histiocytoma Marker Antibody [D11] (orb1252394)

Description

Macrophage / Histiocytoma Marker Antibody [D11]


Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P

Applications

Formalin-fixed,
paraffin-embedded
human ...

Immunogen

A membrane preparation from human hepatocytes was used as the immunogen for the Macrophage / Histiocytoma Marker antibody.

Preservatives

PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide

Form/Appearance


Liquid

Concentration

0.2 mg/mL

Storage

Aliquot and Store at 2-8°C. Avoid freeze-thaw cycles.


Note

For research use only

Application notes

Immunohistochemistry (FFPE): 0.5-1 ug/ml for 30 min at RT (1)Prediluted format: incubate for 30 min at RT (2)Optimal dilution of the Macrophage / Histiocytoma Marker antibody should be determined by the researcher.1. Staining of formalin-fixed tissues requires boiling tissue sections in 10mM Citrate buffer, pH 6.0, for 10-20 min followed by cooling at RT for 20 min.2. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.

Formalin-fixed,
paraffin-embedded
human ...

Isotype

IgG1, kappa

Clonality

Monoclonal

Uniprot ID

NotKnown

Dilution Range

Immunohistochemistry (FFPE): 0.5-1 ug/ml for 30 min at RT (1)Prediluted format: incubate for 30 min at RT (2)Optimal dilution of the Macrophage / Histiocytoma Marker antibody should be determined by the researcher.1. Staining of formalin-fixed tissues requires boiling tissue sections in 10mM Citrate buffer, pH 6.0, for 10-20 min followed by cooling at RT for 20 min.2. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope