

Product Datasheet

ALB Antibody (orb1252418)

Description	ALB Antibody
Species/Host	Mouse
Reactivity	Bovine
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Purified bovine serum albumin was used as the immunogen for the BSA antibody.
Target	ALB
Preservatives	PBS and 0.05% sodium azide
Form/Appearance	Liquid
Concentration	0.2 mg/mL
Storage	Aliquot and Store at 2-8°C. Avoid freez-thaw cycles.
Note	For research use only
Application notes	ELISA: 1-5 ug/ml (order Ab without BSA)Optimal dilution of the BSA antibody should be determined by the researcher.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P02769
Dilution Range	ELISA: 1-5 ug/ml (order Ab without BSA)Optimal dilution of the BSA antibody should be determined by the researcher.
Expiration Date	12 months from date of receipt.