

Product Datasheet

NAPSA Antibody (orb1252434)

Description

NAPSA Antibody

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

FC, IF, IHC-P, WB

Immunogen

A recombinant fragment from the human protein was used as the immunogen for the Napsin A antibody.

Target

NAPSA

Preservatives

PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide

Form/Appearance

Liquid

Concentration

0.2 mg/mL

Storage

Aliquot and Store at 2-8°C. Avoid freeze-thaw cycles.

Note

For research use only

Application notes

Flow Cytometry: 0.5-1 ug/million cells in 0.1mMImmunofluorescence: 1-2 ug/mlWestern blot: 1-2 ug/mlImmunohistochemistry (FFPE): 1-2 ug/ml for 30 min at RT (1)Prediluted format: incubate for 30 min at RT (2)Optimal dilution of the Napsin A antibody should be determined by the researcher.1. Staining of formalin-fixed tissues requires boiling tissue sections in 10mM Citrate buffer, pH 6.0, for 10-20 min followed by cooling at RT for 20 min.2. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.

Isotype

IgG1, kappa

Clonality


Monoclonal

Uniprot ID


O96009

Dilution Range


Flow Cytometry: 0.5-1 ug/million cells in 0.1mMImmunofluorescence: 1-2 ug/mlWestern blot: 1-2 ug/mlImmunohistochemistry (FFPE): 1-2 ug/ml for 30 min at RT (1)Prediluted format: incubate for 30 min at RT (2)Optimal dilution of the Napsin A antibody should be determined by the researcher.1. Staining of formalin-fixed tissues requires boiling tissue sections in 10mM Citrate buffer, pH 6.0, for 10-20 min followed by cooling at RT for 20 min.2. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section



Formalin-fixed, paraffin-embedded human ...



Western blot of 1) K562 and 2) HEK293 ce...



SDS-PAGE Analysis of Purified, BSA-Free ...