

Product Datasheet

HSP90B1 Antibody (orb1252456)

Description

HSP90B1 Antibody

Species/Host

Rat

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IF, IHC-P, WB

Applications
Immunogen

Recombinant full-length human protein was used as the immunogen for the GRP94 antibody.

Target

HSP90B1

Preservatives

PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide

Form/Appearance

Liquid

Concentration

0.2 mg/mL

Storage

Aliquot and Store at 2-8°C. Avoid freeze-thaw cycles.

Note

For research use only

Application notes

Flow Cytometry: 0.5-1 ug/million cells in 0.1ml
Immunofluorescence: 0.5-1 ug/ml
Western blot: 0.5-1 ug/ml
Immunohistology (FFPE): 0.5-1 ug/ml for 30 min at RT
(1) Prediluted format: incubate for 30 min at RT (2) Optimal dilution of the GRP94 antibody should be determined by the researcher.
1. Staining of formalin-fixed tissues requires boiling tissue sections in 10mM Citrate buffer, pH 6.0, for 10-20 min followed by cooling at RT for 20 min.
2. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.

Isotype


IgG2a, kappa

Clonality

Monoclonal

Uniprot ID
P14625
Dilution Range

Flow Cytometry: 0.5-1 ug/million cells in 0.1ml
Immunofluorescence: 0.5-1 ug/ml
Western blot: 0.5-1 ug/ml
Immunohistology (FFPE): 0.5-1 ug/ml for 30 min at RT
(1) Prediluted format: incubate for 30 min at RT (2) Optimal dilution of the GRP94 antibody should be determined by the researcher.
1. Staining of formalin-fixed tissues requires boiling tissue sections in 10mM Citrate buffer, pH 6.0, for 10-20 min followed by cooling at RT for 20 min.
2. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and



Formalin-fixed,
paraffin-embedded
human ...