

Product Datasheet

TLR2 Antibody (orb1252459)

Description

TLR2 Antibody

Species/Host

Mouse

Reactivity

Canine, Human, Monkey

Conjugation

Unconjugated

Tested

FA, FC, IF, IHC-P, WB

Applications

Immunogen

Recombinant human protein was used as the immunogen for the TLR2 antibody.

Target

TLR2

Preservatives

PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide

Form/Appearance

Liquid

Concentration

0.2 mg/mL

Storage

Aliquot and Store at 2-8°C. Avoid freez-thaw cycles.

Note

For research use only

Application notes

Functional Activity (Order SAF formulation)Flow Cytometry: 0.5-1 ug/million cells in 0.1mlImmunofluorescence: 1-2 ug/mlImmunohistochemistry (FFPE): 0.5-1 ug/ml for 30 min at RT (1)Optimal dilution of the TLR2 antibody should be determined by the researcher.1. Staining of formalin-fixed tissues requires boiling tissue sections in 10mM Citrate buffer, pH 6.0, for 10-20 min followed by cooling at RT for 20 min.

Isotype

IgG2a, kappa

Clonality

Monoclonal

Uniprot ID


O60603

Dilution Range

Functional Activity (Order SAF formulation)Flow Cytometry: 0.5-1 ug/million cells in 0.1mlImmunofluorescence: 1-2 ug/mlImmunohistochemistry (FFPE): 0.5-1 ug/ml for 30 min at RT (1)Optimal dilution of the TLR2 antibody should be determined by the researcher.1. Staining of formalin-fixed tissues requires boiling tissue sections in 10mM Citrate buffer, pH 6.0, for 10-20 min followed by cooling at RT for 20 min.

Expiration Date

12 months from date of receipt.



SDS-PAGE
Analysis of
Purified, BSA-
Free ...