

Product Datasheet

TGFA Antibody (orb1252463)

Description

TGFA Antibody

Species/Host

Mouse

Reactivity

Human, Rabbit, Zebrafish

Conjugation

Unconjugated

Tested

FC, IF, IHC-P

Applications

Immunogen

Amino acids 34-43 (PPVAAAVVSH) from the human protein were used as the immunogen for the TGFA antibody.

Target

TGFA

Preservatives

PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide

Form/Appearance

Liquid

Concentration

0.2 mg/mL

Storage

Aliquot and Store at 2-8°C. Avoid freez-thaw cycles.

Note

For research use only

Application notes

Flow Cytometry: 0.5-1 ug/million cells in 0.1ml
Immunofluorescence: 1-2 ug/ml
Immunohistochemistry (FFPE): 2-4 ug/ml for 30 min at RT (1)
Optimal dilution of the TGFA antibody should be determined by the researcher.1.
Staining of formalin-fixed tissues requires boiling tissue sections in 10mM Citrate buffer, pH 6.0, for 10-20 min followed by cooling at RT for 20 min.

Isotype

IgG1, kappa

Clonality

Monoclonal

Uniprot ID


P01135

Dilution Range

Flow Cytometry: 0.5-1 ug/million cells in 0.1ml
Immunofluorescence: 1-2 ug/ml
Immunohistochemistry (FFPE): 2-4 ug/ml for 30 min at RT (1)
Optimal dilution of the TGFA antibody should be determined by the researcher.1.
Staining of formalin-fixed tissues requires boiling tissue sections in 10mM Citrate buffer, pH 6.0, for 10-20 min followed by cooling at RT for 20 min.

Expiration Date

12 months from date of receipt.



Formalin-fixed,
paraffin-embedded
human ...