



Product Datasheet

Slc5a1 antibody (orb125249)

Description	Goat polyclonal antibody to Slc5a1
Species/Host	Goat
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	Slc5a1 (rat)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 70kDa band observed in Rat Kidney lysates (calculated MW of 72.1kDa according to NP_037165.1). Recommended concentration: 0.3-2 µg/ml. Primary incubation 1 hour at room temperature.. IHC: Paraffin embedded Rat Kidney. Recommended concentration: 10 µg/ml.
Protein Sequence	KDTIEIDAEAPQKEK
Clonality	Polyclonal
MW	73.1
Dilution Range	ELISA: 1:32000, WB: 0.3-2 µg/ml, IHC-P: 10 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Rat Kidn...



Immunohistochemical staining of Rat Kidn...



Western blot analysis of Rat Kidney lysa...