


Product Datasheet

ITGB2 Antibody (orb1252505)

Description	ITGB2 Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FA, FC, IF, WB
Immunogen	Stimulated human leukocytes were used as the immunogen for the CD18 antibody.
Target	ITGB2
Preservatives	PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide
Form/Appearance	Liquid
Concentration	0.2 mg/mL
Storage	Aliquot and Store at 2-8°C. Avoid freez-thaw cycles.
Note	For research use only
Application notes	Functional Activity (Order SAF formulation)Flow Cytometry: 0.5-1 ug/million cells in 0.1mlImmunofluorescence: 0.5-1 ug/mlOptimal dilution of the CD18 antibody should be determined by the researcher.
Isotype	IgG2a, kappa
Clonality	Monoclonal
Uniprot ID	P05107
Dilution Range	Functional Activity (Order SAF formulation)Flow Cytometry: 0.5-1 ug/million cells in 0.1mlImmunofluorescence: 0.5-1 ug/mlOptimal dilution of the CD18 antibody should be determined by the researcher.
Expiration Date	12 months from date of receipt.



SDS-PAGE
Analysis of
Purified, BSA-
Free ...