


Product Datasheet

ICAM1 Antibody (orb1252523)

Description	ICAM1 Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FA, FC, IF, WB
Immunogen	Human umbilical cord vein endothelial cells (HUVEC) were used as the immunogen for the ICAM1 antibody.
Target	ICAM1
Preservatives	PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide
Form/Appearance	Liquid
Concentration	0.2 mg/mL
Storage	Aliquot and Store at 2-8°C. Avoid freez-thaw cycles.
Note	For research use only
Application notes	Flow Cytometry: 0.5-1 ug/million cells in 0.1ml Immunofluorescence: 0.5-1 ug/ml Functional Activity (Order SAF formulation) Optimal dilution of the ICAM1 antibody should be determined by the researcher.
Isotype	IgG2a, kappa
Clonality	Monoclonal
Uniprot ID	P05362
Dilution Range	Flow Cytometry: 0.5-1 ug/million cells in 0.1ml Immunofluorescence: 0.5-1 ug/ml Functional Activity (Order SAF formulation) Optimal dilution of the ICAM1 antibody should be determined by the researcher.
Expiration Date	12 months from date of receipt.



SDS-PAGE
Analysis of
Purified, BSA-
Free ...