


Product Datasheet

HLA-DRB1 Antibody (orb1252528)

Description	HLA-DRB1 Antibody
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Primate, Rodent
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Activated human peripheral blood mononuclear cells were used as the immunogen for the HLA-DRB1 antibody.
Target	HLA-DRB1
Preservatives	PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide
Form/Appearance	Liquid
Concentration	0.2 mg/mL
Storage	Aliquot and Store at 2-8°C. Avoid freez-thaw cycles.
Note	For research use only
Application notes	Flow Cytometry: 0.5-1 ug/million cells in 0.1ml Immunofluorescence: 0.5-1 ug/ml Optimal dilution of the HLA-DRB1 antibody should be determined by the researcher.
Isotype	IgG2a, kappa
Clonality	Monoclonal
Uniprot ID	P01911
Dilution Range	Flow Cytometry: 0.5-1 ug/million cells in 0.1ml Immunofluorescence: 0.5-1 ug/ml Optimal dilution of the HLA-DRB1 antibody should be determined by the researcher.
Expiration Date	12 months from date of receipt.



SDS-PAGE
Analysis of
Purified, BSA-
Free ...