


Product Datasheet

HLA-ABC Antibody [246-B8.E7] (orb1252533)

Description	HLA-ABC Antibody [246-B8.E7]
Species/Host	Mouse
Reactivity	Human, Mouse, Primate
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Human peripheral blood leukocytes were used as the immunogen for the HLA-ABC antibody.
Preservatives	PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide
Form/Appearance	Liquid
Concentration	0.2 mg/mL
Storage	Aliquot and Store at 2-8°C. Avoid freez-thaw cycles.
Note	For research use only
Application notes	Western blot: 1-2 ug/mlFlow Cytometry: 0.5-1 ug/million cells in 0.1mlImmunofluorescence: 0.5-1 ug/mlOptimal dilution of the HLA-ABC antibody should be determined by the researcher.
Isotype	IgG2a, kappa
Clonality	Monoclonal
Uniprot ID	P30443, P30499, P01889
Dilution Range	Western blot: 1-2 ug/mlFlow Cytometry: 0.5-1 ug/million cells in 0.1mlImmunofluorescence: 0.5-1 ug/mlOptimal dilution of the HLA-ABC antibody should be determined by the researcher.
Expiration Date	12 months from date of receipt.



Western blot testing of ThP-1 cell lysat...



SDS-PAGE Analysis of Purified, BSA-Free ...