

Product Datasheet

PCBP3 antibody (orb125257)

Description	Goat Polyclonal to PCBP3. NP_065389.2; NP_001123613.1, PMID: 14612387 and GenelD: 54039.
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Target	PCBP3 (aa14-28)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Preliminary experiments in Human Testis, Adipose, Lung and several cel line lysates gave no specific signal but low background (at antibody concentration up to 2µg/ml). We would appreciate any feedback from people in the field - have any results been repo
Protein Sequence	PHSTLSTLSHHPQPQ
Clonality	Polyclonal
MW	39.5; 36.8
Entrez	54039
Dilution Range	ELISA: 1:32000, WB: 2 µg/ml
Expiration Date	12 months from date of receipt.