

Product Datasheet

CNN1 Antibody (orb1252578)

Description

CNN1 Antibody

Species/Host

Mouse

Reactivity

Human, Rat

Conjugation

Unconjugated

Tested Applications

FC, IF, IHC-P, WB

Immunogen

Crude human uterus extract was used as the immunogen for the Calponin antibody.

Target

CNN1

Preservatives

PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide

Form/Appearance

Liquid

Concentration

0.2 mg/mL

Storage

Aliquot and Store at 2-8°C. Avoid freez-thaw cycles.

Note

For research use only

Application notes

Flow Cytometry: 0.5-1 ug/million cells in 0.1mMImmunofluorescence: 1-2 ug/mlImmunohistochemistry (FFPE): 0.5-1 ug/ml for 30 min at RT (1)Optimal dilution of the Calponin antibody should be determined by the researcher.1. Staining of formalin-fixed tissues requires boiling tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 10-20 min followed by cooling at RT for 20 min

Isotype

IgG1, kappa

Clonality

Monoclonal

Uniprot ID


P51911

Dilution Range


Flow Cytometry: 0.5-1 ug/million cells in 0.1mMImmunofluorescence: 1-2 ug/mlImmunohistochemistry (FFPE): 0.5-1 ug/ml for 30 min at RT (1)Optimal dilution of the Calponin antibody should be determined by the researcher.1. Staining of formalin-fixed tissues requires boiling tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 10-20 min followed by cooling at RT for 20 min

Expiration Date

12 months from date of receipt.



Formalin-fixed, paraffin-embedded human ...



SDS-PAGE Analysis of Purified, BSA-Free ...