

Product Datasheet

CHGA Antibody (orb1252582)

Description
CHGA Antibody
Species/Host

Mouse

Reactivity

Human, Rat

Conjugation

Unconjugated

Tested

FC, IF, IHC-P

Applications
Immunogen

Recombinant human protein was used as the immunogen for the Chromogranin A antibody.

Target

CHGA

Preservatives

PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide

Form/Appearance

Liquid

Concentration

0.2 mg/mL

Storage

Aliquot and Store at 2-8°C. Avoid freeze-thaw cycles.

Note

For research use only

Application notes

Flow Cytometry: 0.5-1 ug/million cells in 0.1ml
Immunofluorescence: 0.5-1 ug/ml
Immunohistochemistry (FFPE): 0.25-0.5 ug/ml for 30 min at RT (1)Prediluted format : incubate for 30 min at RT (2)Optimal dilution of the Chromogranin A antibody should be determined by the researcher.1. Staining of formalin-fixed tissues requires boiling tissue sections in 10mM Citrate buffer, pH 6.0, for 10-20 min followed by cooling at RT for 20 minutes2. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.

Isotype


IgG1, kappa

Clonality


Monoclonal

Uniprot ID
P10645
Dilution Range


Flow Cytometry: 0.5-1 ug/million cells in 0.1ml
Immunofluorescence: 0.5-1 ug/ml
Immunohistochemistry (FFPE): 0.25-0.5 ug/ml for 30 min at RT (1)Prediluted format : incubate for 30 min at RT (2)Optimal dilution of the Chromogranin A antibody should be determined by the researcher.1. Staining of formalin-fixed tissues requires boiling tissue sections in 10mM Citrate buffer, pH 6.0, for 10-20 min followed by cooling at RT for 20 minutes2. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section



Formalin-fixed, paraffin-embedded human ...



Formalin-paraffin human pancreas stained...



Formalin-paraffin rat pancreas stained w...