


Product Datasheet

SEMA4D Antibody (orb1252586)

Description	SEMA4D Antibody
Species/Host	Mouse
Reactivity	Human, Monkey, Mouse
Conjugation	Unconjugated
Tested Applications	FA, FC, IF, WB
Immunogen	PHA stimulated human peripheral blood lymphocytes were used as the immunogen for the Semaphorin-4D antibody.
Target	SEMA4D
Preservatives	PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide
Form/Appearance	Liquid
Concentration	0.2 mg/mL
Storage	Aliquot and Store at 2-8°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	Functional Activity (Order SAF formulation)Flow Cytometry: 0.5-1 ug/million cells in 0.1mlImmunofluorescence: 0.5-1 ug/mlOptimal dilution of the Semaphorin-4D antibody to be determined by the researcher.
Isotype	IgG1, kappa
Clonality	Monoclonal
Uniprot ID	Q92854
Dilution Range	Functional Activity (Order SAF formulation)Flow Cytometry: 0.5-1 ug/million cells in 0.1mlImmunofluorescence: 0.5-1 ug/mlOptimal dilution of the Semaphorin-4D antibody to be determined by the researcher.
Expiration Date	12 months from date of receipt.



SDS-PAGE
Analysis of
Purified, BSA-
Free ...