

Product Datasheet

MS4A1 Antibody (orb1252625)

Description

MS4A1 Antibody

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IF, IHC

Applications

Immunogen

Recombinant human MS4A1/CD20 protein was used as the immunogen for this anti-CD20 antibody.

Target

MS4A1

Preservatives

PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide

Form/Appearance

Liquid

Concentration

0.2 mg/mL

Storage

Aliquot and Store at 2-8°C. Avoid freez-thaw cycles.

Note

For research use only

Application notes

Flow Cytometry: 0.5-1 ug/million cells
IF: 0.5-1 ug/ml
Immunohistology (FFPE): 0.5-1 ug/ml for 30 minutes at RT
(1) Prediluted format : incubate for 30 min at RT
(2) The concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the anti-CD20 antibody to be titrated up or down for optimal performance.
1. Staining of formalin-fixed tissues requires boiling tissue sections in 10mM citrate buffer, pH 6.0, for 10-20 min followed by cooling at RT for 20 minutes.
2. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.

Isotype

IgG2a, kappa

Clonality


Monoclonal

Uniprot ID

P11836

Dilution Range

Flow Cytometry: 0.5-1 ug/million cells
IF: 0.5-1 ug/ml
Immunohistology (FFPE): 0.5-1 ug/ml for 30 minutes at RT
(1) Prediluted format : incubate for 30 min at RT
(2) The concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the anti-CD20 antibody to be titrated up or down for optimal performance.
1. Staining of formalin-fixed tissues requires boiling tissue sections in 10mM citrate buffer, pH 6.0, for 10-20 min followed by



IHC testing of anti-CD20 antibody and FFPE...

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0) 1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina 27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Expiration Date

12 months from date of receipt.