


Product Datasheet

CD4 Antibody (orb1252634)

Description	CD4 Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Recombinant human CD4 protein was used as immunogen for the CD4 antibody.
Target	CD4
Preservatives	PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide
Form/Appearance	Liquid
Concentration	0.2 mg/mL
Storage	Aliquot and Store at 2-8°C. Avoid freez-thaw cycles.
Note	For research use only
Application notes	Flow Cytometry: 0.5-1ug/10e6 cellsImmunofluorescence: 0.5-1 ug/mlOptimal dilution of the CD4 antibody for a specific application should be determined by researcher.
Isotype	IgG1, kappa
Clonality	Monoclonal
Uniprot ID	P01730
Dilution Range	Flow Cytometry: 0.5-1ug/10e6 cellsImmunofluorescence: 0.5-1 ug/mlOptimal dilution of the CD4 antibody for a specific application should be determined by researcher.
Expiration Date	12 months from date of receipt.



SDS-PAGE
Analysis of
Purified, BSA-
Free ...