

Product Datasheet

Her2 Antibody (orb1252637)

Description

Her2 Antibody

Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF
Immunogen	Recombinant human HER2 protein was used as the immunogen for this antibody.(1)
Target	Her2
Preservatives	PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide
Form/Appearance	Liquid
Concentration	0.2 mg/mL
Storage	Aliquot and Store at 2-8°C. Avoid freez-thaw cycles.
Note	For research use only
Application notes	Flow Cytometry: 0.5-1 ug/million cellsIF: 0.5-1 ug/mlThe concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the antibody to be titrated up or down for optimal performance.1. This HER2 antibody binds to the extracellular/cell surface region of the protein.
Isotype	IgG1, kappa
Clonality	Monoclonal
Uniprot ID	P04626
Dilution Range	Flow Cytometry: 0.5-1 ug/million cellsIF: 0.5-1 ug/mlThe concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the antibody to be titrated up or down for optimal performance.1. This HER2 antibody binds to the extracellular/cell surface region of the protein.
Expiration Date	12 months from date of receipt.