

Product Datasheet

Nuclear Membrane Marker Antibody [NM97] (orb1252652)

Description

Nuclear Membrane Marker Antibody [NM97]

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IF, WB

Applications

Immunogen

Nuclei of myeloid leukemia biopsy cells were used as the immunogen for this Nuclear membrane marker antibody.

Preservatives

PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide

Form/Appearance

Liquid

Concentration

0.2 mg/mL

Storage

Aliquot and Store at 2-8°C. Avoid freeze-thaw cycles.

Note

For research use only

Application notes

IF: 0.5-1.0 ug/ml
Flow Cytometry: 0.5-1 ug/million cells in 0.1ml
The concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the Nuclear membrane marker antibody to be titrated up or down for optimal performance.

Isotype

IgG1, kappa

Clonality


Monoclonal

Dilution Range

IF: 0.5-1.0 ug/ml
Flow Cytometry: 0.5-1 ug/million cells in 0.1ml
The concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the Nuclear membrane marker antibody to be titrated up or down for optimal performance.

Expiration Date

12 months from date of receipt.



SDS-PAGE
Analysis of
Purified,
BSA-Free ...