

Product Datasheet

Helicobacter Pylori Antibody (orb1252660)

Description

Helicobacter Pylori Antibody

Species/Host

Rabbit

Reactivity

H. pylori

Conjugation

Unconjugated

Tested

FC, IF, IHC

Applications
Immunogen

Total sonicate of Helicobacter pylori was used as the immunogen for this antibody.

Preservatives

PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide

Form/Appearance

Liquid

Concentration

0.2 mg/mL

Storage

Aliquot and Store at 2-8°C. Avoid freeze-thaw cycles.

Note

For research use only

Application notes

Flow Cytometry: 0.5-1 ug/million cells
IF: 1-2 ug/ml
IHC (FFPE): 1-2 ug/ml for 30 minutes at RT
(1)The concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the antibody to be titrated up or down for optimal performance.
1. Staining of formalin-fixed tissues requires boiling tissue sections in 10mM Citrate, pH 6.0, for 10-20 min followed by cooling at RT for 20 minutes.

Isotype

IgG

Clonality


Polyclonal

Dilution Range

Flow Cytometry: 0.5-1 ug/million cells
IF: 1-2 ug/ml
IHC (FFPE): 1-2 ug/ml for 30 minutes at RT
(1)The concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the antibody to be titrated up or down for optimal performance.
1. Staining of formalin-fixed tissues requires boiling tissue sections in 10mM Citrate, pH 6.0, for 10-20 min followed by cooling at RT for 20 minutes.

Expiration Date

12 months from date of receipt.



Formalin/paraffin human stomach stained ...