


Product Datasheet

Cytomegalovirus Antibody [CMV100] (orb1252666)

Description	Cytomegalovirus Antibody [CMV100]
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, WB
Preservatives	PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide
Form/Appearance	Liquid
Concentration	0.2 mg/mL
Storage	Aliquot and Store at 2-8°C. Avoid freez-thaw cycles.
Note	For research use only
Application notes	ELISA: order BSA-free CMV antibody format for coatingImmunocytochemistry (Acetone or paraformaldehyde fixed) 0.5-1 ug/ml for 30 minIF: 0.5-1 ug/mlThe concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the CMV antibody to be titrated up or down for optimal performance.
Isotype	IgG1, kappa
Clonality	Monoclonal
Dilution Range	ELISA: order BSA-free CMV antibody format for coatingImmunocytochemistry (Acetone or paraformaldehyde fixed) 0.5-1 ug/ml for 30 minIF: 0.5-1 ug/mlThe concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the CMV antibody to be titrated up or down for optimal performance.
Expiration Date	12 months from date of receipt.



SDS-PAGE
Analysis of
Purified,
BSA-Free ...