

Product Datasheet

ITGAL Antibody (orb1252745)

Description

ITGAL Antibody

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

FA, FC, IF, WB

Applications

Immunogen

Stimulated human leukocytes were used as the immunogen for this CD11b antibody.

Target

ITGAL

Preservatives

PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide

Form/Appearance

Liquid

Concentration

0.2 mg/mL

Storage

Aliquot and Store at 2-8°C. Avoid freeze-thaw cycles.

Note

For research use only

Application notes

Functional Testing: order BSA/azide-free formatFlow Cytometry: 0.5-1 ug/million cellsIF: 0.5-1 ug/mlThe concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the antibody to be titrated up or down for optimal performance.

Isotype

IgG2b, kappa

Clonality

Monoclonal

Uniprot ID


[P20701](#)

Dilution Range


Functional Testing: order BSA/azide-free formatFlow Cytometry: 0.5-1 ug/million cellsIF: 0.5-1 ug/mlThe concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the antibody to be titrated up or down for optimal performance.

Expiration Date

12 months from date of receipt.



FACS staining of human PBMCs using CD11a...



SDS-PAGE Analysis of Purified, BSA-Free ...