

Product Datasheet

IGL Antibody (orb1252749)

Description

IGL Antibody

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC, WB

Applications

Immunogen

Purified human Ig was used as the immunogen for this Lambda light chain antibody.

Target

IGL

Preservatives

PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide

Form/Appearance

Liquid

Concentration

0.2 mg/mL

Storage

Aliquot and Store at 2-8°C. Avoid freeze-thaw cycles.

Note

For research use only

Application notes

Flow Cytometry: 0.5-1 ug/million cells
IHC (FFPE): 0.5-1 ug/ml for 30 minutes at RT (1)
The concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the antibody to be titrated up or down for optimal performance.
1. Staining of formalin-fixed tissues requires boiling tissue sections in 10mM Citrate Buffer, pH 6.0, for 10-20 min followed by cooling at RT for 20 minutes.

Isotype

IgG1, kappa

Clonality

Monoclonal

Uniprot ID


[P01842](#), [P01701](#)

Dilution Range


Flow Cytometry: 0.5-1 ug/million cells
IHC (FFPE): 0.5-1 ug/ml for 30 minutes at RT (1)
The concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the antibody to be titrated up or down for optimal performance.
1. Staining of formalin-fixed tissues requires boiling tissue sections in 10mM Citrate Buffer, pH 6.0, for 10-20 min followed by cooling at RT for 20 minutes.

Expiration Date

12 months from date of receipt.



Formalin/paraffin human tonsil stained w...



SDS-PAGE Analysis of Purified, BSA-Free ...