


Product Datasheet

CCNA2 Antibody (orb1252820)

Description	CCNA2 Antibody
Species/Host	Mouse
Reactivity	Bovine, Human, Mink, Monkey, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length bovine protein was used as the immunogen for this CCNA2 antibody.
Target	CCNA2
Preservatives	PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide
Form/Appearance	Liquid
Concentration	0.2 mg/mL
Storage	Aliquot and Store at 2-8°C. Avoid freez-thaw cycles.
Note	For research use only
Application notes	Flow Cytometry: 0.5-1 ug/million cellsIF: 1-2 ug/mlThe concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the Cyclin A2 antibody to be titrated up or down for optimal performance.
Isotype	IgG2a, kappa
Clonality	Monoclonal
Uniprot ID	P20248
Dilution Range	Flow Cytometry: 0.5-1 ug/million cellsIF: 1-2 ug/mlThe concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the Cyclin A2 antibody to be titrated up or down for optimal performance.
Expiration Date	12 months from date of receipt.



SDS-PAGE
Analysis of
Purified, BSA-
Free ...