


Product Datasheet

C5ar2 Antibody (orb1252848)

Description	C5ar2 Antibody
Species/Host	Rabbit
Reactivity	Human, Primate
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A recombinant protein fragment from the N-terminal region of human GPR77 was used as the immunogen for this antibody.
Target	C5ar2
Preservatives	PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Aliquot and Store at -20°C. Avoid freez-thaw cycles.
Note	For research use only
Application notes	Western blot: 1-3 ug/ml Provided assay concentrations are suggestions only, GPR77 antibody titration may be required for optimal results.
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q695P6
Dilution Range	Western blot: 1-3 ug/ml Provided assay concentrations are suggestions only, GPR77 antibody titration may be required for optimal results.
Expiration Date	12 months from date of receipt.



Western blot testing of human samples us...