



Product Datasheet

IPO9 Antibody (orb1252855)

Description	IPO9 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A recombinant protein fragment from the C-terminal region of human Importin-9 was used as the immunogen for this antibody.
Target	IPO9
Preservatives	PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Aliquot and Store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	Western blot: 0.25-1 ug/ml Provided assay concentrations are suggestions only, Importin-9 antibody titration may be required for optimal results.
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96P70
Dilution Range	Western blot: 0.25-1 ug/ml Provided assay concentrations are suggestions only, Importin-9 antibody titration may be required for optimal results.
Expiration Date	12 months from date of receipt.



Western blot testing of human samples wi...



Western blot testing of mouse samples wi...