


Product Datasheet

NTSR2 Antibody (orb1252867)

Description	NTSR2 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A recombinant protein fragment from the N-terminal half of human Neurotensin Receptor 2 was used as the immunogen for this antibody.
Target	NTSR2
Preservatives	PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Aliquot and Store at -20°C. Avoid freez-thaw cycles.
Note	For research use only
Application notes	Western blot: 1-3 ug/ml Provided assay concentrations are suggestions only, Neurotensin Receptor 2 antibody titration may be required for optimal results.
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O95665
Dilution Range	Western blot: 1-3 ug/ml Provided assay concentrations are suggestions only, Neurotensin Receptor 2 antibody titration may be required for optimal results.
Expiration Date	12 months from date of receipt.



Western blot testing of human samples wi...