

Product Datasheet

TDRD9 antibody (orb1253)

Description	TDRD9 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr
Immunogen	KLH conjugated synthetic peptide derived from human TDRD9/HIG1 (801-900/1382 aa)
Target	TDRD9
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	155 kDa
Uniprot ID	Q14BI7
Expiration Date	12 months from date of receipt.