

Product Datasheet

PXN Antibody (orb1253069)

Description	PXN Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Paxillin (phospho Ser272) antibody was raised against a peptide sequence around phosphorylation site of Serine 272 (M-A-S (p)-L-S) derived from Human Paxillin.
Target	PXN
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Application notes	Western Blot: 1:500~1:1000
Isotype	IgG
Clonality	Polyclonal
MW	68 kDa
Uniprot ID	P49023
NCBI	NP_001074324.1
Dilution Range	Western Blot: 1:500~1:1000
Expiration Date	12 months from date of receipt.