


Product Datasheet

STK39 Antibody (orb1253076)

Description	STK39 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Peptide sequence around phosphorylation site of Serine 309(K-G-S(p)-F-R) derived from Human STK39.
Target	STK39
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Application notes	Western Blot: 1:500~1:1000, Immunohistochemistry: 1:50~1:100
Isotype	IgG
Clonality	Polyclonal
MW	60 kDa
Uniprot ID	Q9UEW8
NCBI	NP_037365.2
Dilution Range	Western Blot: 1:500~1:1000, Immunohistochemistry: 1:50~1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...