


Product Datasheet

EIF2AK3 Antibody (orb1253080)

Description	EIF2AK3 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	PEK/PERK (phospho Thr982) antibody was raised against a peptide sequence around phosphorylation site of threonine 981 (A-R-H (p)-T-G) derived from Human PEK/PERK.
Target	EIF2AK3
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	125 kDa
Uniprot ID	Q9NZJ5
NCBI	NP_004827.4
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...