


Product Datasheet

OXSR1 Antibody (orb1253083)

Description

OXSRI Antibody



Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

IHC, WB

Applications

Immunohistochemical
analysis of paraffin...

Immunogen

OSR1 (phospho Thr185) antibody was raised against a peptide sequence around phosphorylation site of threonine185 (R-K-T (p)-F-V) derived from Human OSR1.

Target

OXSRI

Preservatives

Antibody supplied in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

Store antibody at -20°C for up to one year.

Note

For research use only

Application notes

Western Blot: 1:500~1:1000, Immunohistochemistry: 1:50~1:100

Isotype

IgG

Clonality

Polyclonal

MW

65 kDa

Uniprot ID

[O95747](#)

NCBI

[NP_005100.1](#)

Dilution Range

Western Blot: 1:500~1:1000, Immunohistochemistry: 1:50~1:100

Expiration Date

12 months from date of receipt.