


Product Datasheet

MERTK, TYRO3 Antibody (orb1253089)

Description

MERTK, TYRO3 Antibody



Species/Host

Rabbit

Reactivity

Human, Mouse, Rat


Conjugation

Unconjugated

Tested

IHC, WB

Applications



Immunohistochemical
analysis of paraffin...

Immunogen

MER/SKY (phospho Tyr749/681) antibody was raised against a peptide sequence around phosphorylation site of tyrosine 749/681 (K-I-Y (p)-S-G) derived from Human MER/SKY .

Target

MERTK, TYRO3

Preservatives

Antibody supplied in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

Store antibody at -20°C for up to one year.

Note

For research use only

Application notes

Western Blot: 1:500~1:1000,
Immunohistochemistry: 1:50~1:100

Isotype

IgG

Clonality

Polyclonal

MW

97 kDa

Uniprot ID

[Q06418](#), [Q12866](#)

NCBI

[NP_006334.2](#)

Dilution Range

Western Blot: 1:500~1:1000,
Immunohistochemistry: 1:50~1:100

Expiration Date

12 months from date of receipt.