

Product Datasheet

TAU Antibody (orb1253104)

Description	TAU Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Tau (phospho Thr534/217) antibody was raised against a peptide sequence around phosphorylation site of threonine 534/217 (L-P-T (p)-P-P) derived from Human Tau.
Target	TAU
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Application notes	Western Blot: 1:500~1:1000
Isotype	IgG
Clonality	Polyclonal
MW	65 kDa
Uniprot ID	P10636
NCBI	NP_058519.3
Dilution Range	Western Blot: 1:500~1:1000
Expiration Date	12 months from date of receipt.