


Product Datasheet

CHEK2 Antibody (orb1253137)

Description	CHEK2 Antibody	
Species/Host	Rabbit	
Reactivity	Human, Mouse, Rat	
Conjugation	Unconjugated	
Tested Applications	IF, WB	Immunofluorescence staining of methanol-...
Immunogen	Chk2 (phospho Thr383) antibody was raised against a peptide sequence around phosphorylation site of threonine383 (M-R-T (p)-L-C) derived from Human Chk1.	
Target	CHEK2	
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.	
Form/Appearance	Liquid	
Concentration	1 mg/mL	
Storage	Store antibody at -20°C for up to one year.	
Note	For research use only	
Application notes	Western Blot: 1:500~1:1000, Immunofluorescence: 1:100~1:200	
Isotype	IgG	
Clonality	Polyclonal	
MW	60 kDa	
Uniprot ID	O96017	
NCBI	NP_009125.1	
Dilution Range	Western Blot: 1:500~1:1000, Immunofluorescence: 1:100~1:200	
Expiration Date	12 months from date of receipt.	