



Product Datasheet

PMS2 antibody (orb125320)

Description	Rabbit Polyclonal to PMS2.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P54278, Q68D20
Dilution Range	WB:1: 500-1: 3000, IHC-P: 1:50-1:100, ELISA: 1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from 3...



Immunohistochemical analysis of formalin...