


Product Datasheet

FOXP2 Antibody (orb1253230)

Description	FOXP2 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	FOXP2 antibody was raised against a peptide sequence around aa.707~711 (E-E-P-L-S) derived from Human FOXP2.
Target	FOXP2
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Application notes	Western Blot: 1:500~1:1000
Clonality	Polyclonal
MW	85 kDa
Uniprot ID	O15409
NCBI	NP_055306.1
Dilution Range	Western Blot: 1:500~1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of extract from JK...