


Product Datasheet

GRIA2 Antibody (orb1253333)

Description	GRIA2 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Glutamate receptor 2 (Precursor) (Ab-880) antibody was raised against a peptide sequence around aa.878~882 (I-E-S-V-K) derived from Human Glutamate receptor 2.
Target	GRIA2
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Application notes	Western Blot: 1:500~1:1000
Clonality	Polyclonal
MW	100 kDa
Uniprot ID	P42262
NCBI	NP_000817.2
Dilution Range	Western Blot: 1:500~1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of lysed extracts ...