



Product Datasheet

RAB3GAP1 antibody (orb125334)


Description	Rabbit Polyclonal to RAB3GAP1.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q15042
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, ELISA: 1:40000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from 2...



Western blot analysis of extracts from J...



Immunohistochemical analysis of formalin...