



Product Datasheet

PTK2 Antibody (orb1253429)

Description	PTK2 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	FAK (phospho-Tyr576/Tyr577) antibody was raised against a peptide sequence around phosphorylation site of tyrosine 576/tyrosine 577 (S-T-Y (p) -Y (p) -K-A) derived from Human FAK.
Target	PTK2
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Application notes	Western Blot: 1:500~1:1000, Immunofluorescence: 1:100~1:200
Clonality	Polyclonal
MW	125 kDa
Uniprot ID	Q05397
NCBI	NP_005598.3
Dilution Range	Western Blot: 1:500~1:1000, Immunofluorescence: 1:100~1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of lysed extracts ...



Immunofluorescence staining of methanol-...