


Product Datasheet

CTTN Antibody (orb1253451)

Description	CTTN Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Cortactin (Phospho-Tyr466) antibody was raised against a peptide sequence around phosphorylation site of tyrosine 466 (P-V-Y (p) -E-T) derived from Human CORTACTIN.
Target	CTTN
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Application notes	Western Blot: 1:500~1:1000
Clonality	Polyclonal
MW	85 kDa
Uniprot ID	Q14247
NCBI	NP_005222.2
Dilution Range	Western Blot: 1:500~1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of lysed extracts ...